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
User guide

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Overview

Vaginal weights strengthen the pelvic floor muscles discreetly, without requiring additional time for exercises. The set includes four beads of three different weights for gradual progression.

Instructions for Use

To start, wear the vaginal weights in a place where you can easily remove or adjust them. Begin by wearing them for about 30 minutes a day to get used to their presence. Your position and activity (such as sitting or walking) influence the difficulty. Start by sitting, then gradually progress to more demanding activities like household chores or squats.

Use light weights at first and gradually increase the difficulty based on your comfort level. Consistency is key; there is no standard duration for their use.

Insertion

Choose a comfortable position: squatting, sitting, or lying down. Insert the weights until the vaginal entrance (about 2 cm deep) closes, indicating proper placement. If inserted too deeply, the pelvic muscles may not hold them effectively. To facilitate insertion, use a water-based vaginal lubricant and ensure the retrieval strap remains outside the vagina for easy removal.

Removal

Gently pull on the retrieval strap to remove the beads from the vagina. For beads without a strap, relax your pelvic muscles to allow them to come out naturally. If removal is difficult, relax and apply a water-based lubricant, then try standing up, walking, coughing, or jumping in place. You can also insert a finger behind the bead(s) inside the vagina and gently pull downward.


Maintenance

Remove the beads from the strap to clean the vaginal weights before and after each use with a cleaner such as Steri-clean. Store the weights in their original packaging.

Warranty

For more information about your product's warranty, please visit our website: floravi.com.

For any questions or feedback, contact us at: info@floravi.com

Single weights / Poids seuls

30g 40g 50g

Multiple weights / Poids combinés

30g + 30g = 60g 30g + 40g = 70g
30g + 50g = 80g 40g + 50g = 90g