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
User guide

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## Overview

Vaginal dilators help relax the pelvic floor muscles and can be useful if you experience difficulties during sexual intercourse, tampon insertion, or gynecological exams. They are designed to improve comfort and contribute to a more satisfying sexual life.

## Insertion

For optimal comfort, warm the dilator with hot water, relax, and apply a water-based lubricant to ease insertion. Choose a comfortable position: lying on your back, sitting with your knees to your chest, or standing with your knees slightly flexed.

Gently insert the appropriate dilator for your needs (it's better to underestimate the size). Allow your vagina to get accustomed to its presence before introducing movement. Slowly incorporate circular, side-to-side, and in-and-out motions.

If you experience pain during removal, stop, take a deep breath, and cough gently. When you're ready to move to a larger size, start with the previous dilator—the one that caused no discomfort—to facilitate the transition.

## Usage

The results of using dilators vary from person to person; regular practice helps achieve benefits faster. However, it is essential to listen to your body and progress at your own pace. A slight tingling sensation when transitioning to a larger size is normal and should decrease over time. If the sensation persists, return to a smaller dilator.

## Maintenance

For hygiene reasons, clean the dilator before and after each use with a cleaner, such as Steri-clean. Store it in its original packaging to extend its lifespan. Sharing dilators is strongly discouraged.

## Warnings

Use water-based lubricant. If you notice any bleeding, discontinue use and consult a healthcare professional.

## Dimensions